FORM	NO	



CENTRALRESEARCHFACILITY INDIANINSTITUTEOFTECHNOLOGYKHARAGPUR

REQUISITION FORM FORUSE OF AFM

I would like to use the following instrument facility for analysis of my samples. I do also declare that I shall acknowledge the relevant facility in CRF in the future Publications.

SampleSizeUseseparateformsfordifferentanalysis/instrumentsunderdownloadlink		
Name of the Supervisor:		
Name of the User's:		
EC/Registration/Roll:		
Department /Centre /School:		
Contact No. :		
E-mail ID:		

Instructions for Sample Preparation:

- 1. No loose powders, corrosive samples, liquids/vapours will be done.
- 2. One sample will roughly take 1 hour for imaging. Maximum three samples (same material) will be done in a 3 hours slot. Only two samples will be done if samples are of different materials.
- 3. Please DO NOT bring samples wrapped in tissue paper/ foil. Please bring samples in petry dish. Bring your own tweezers. Do not touch sample with hand.
- 4. Normally imaging in intermittent contact mode (tapping mode) will be performed. You might have to wait longer for contact mode imaging.
- 5. Please do not bring very rough (roughness ~ 10 micron or higher) samples.
- 6. Maximum sample thickness is 5 mm.
- 7. Not all modes are available at any given time due to either set up constraints or other reasons.
- 8.For AFM max sample size is 15 mm x 15 mm and max sample thickness is < 5 mm.

TERMS&CONDITIONSANDCHARGESFORAVAILINGINSTRUMENTFACILITIES

- $1. \quad For the times lot and other necessary criteria, please contact the concerned staff beforehand$
- 2. Repeat of analysis may be done if there is any instrumental or operating error.
- 3. Acknowledge the CRF Facility in the future publication if any.

CHARGES	DURATION	AMOUNTININR
Internal(Institute)	1slot=3Hrs	750/-
Internal(project)	1slot=3Hrs	900/-
External users(academic)	1slot=3Hrs	3000+ 18%GST
R&D labs /Industry	1slot=3Hrs	5000+ 18% GST